LANGE 

DOC022.98.90373

LZV936
DR 6000

02/2012, Edition 1A

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General information

Safety notes

Read through the entire user manual carefully before you unpack the device, set up and put into operation. Observe all danger and warning notes. Non-observance could lead to serious injury of the operator or to damage to the device.

To ensure safety when using the instrument, the DR6000 may not be used in any manner that is not described in this user manual.

Danger notes in this user manual

 DANGER
Indicates a potentially or imminently dangerous situation that, if it is not avoided, leads to death or to serious injuries.

 WARNING
Indicates a potentially or imminently dangerous situation that, if it is not avoided, can lead to death or to serious injuries.

 CAUTION
Indicates a possible dangerous situation that can have minor or moderate injuries as the result.

NOTICE
Indicates a situation that, if it is not avoided, can lead to damage to the device. Information that particularly should be emphasized.

Note: Information that supplements aspects from the main text.

Warning labels

Observe all marks and labels that are attached to the device. Non-observance may result in personal injury or damage to the device. For symbols attached to the instrument, corresponding warning notes are found in the user manual.

	This symbol may be attached to the device and references the operation- and/or safety notes in the user manual.
	Electrical equipment marked with this symbol may as of August 12, 2005 Europe-wide no longer be disposed of in unsorted house or industrial waste. According to valid provisions (EU Directive 2002/96/EC), from this point consumers in the EU must return old electrical devices to the manufacturer for disposal. This is free for the consumer. Note: Contact the manufacturer or supplier to find out how you can return worn out devices, electrical accessories supplied by the manufacturer and all auxiliary articles for correct disposal or recycling.

Installation notes

1. Select **SYSTEM CHECKS > INSTRUMENT UPDATE**.
2. Connect the USB stick with the additional software "Brewery Analysis LZV936" to a USB connection (type A) of the DR 6000.
3. Confirm with **OK**.
The software is installed.
4. Switch the DR 6000 off after the complete installation at the power source and start it again.
Note: Wait at least **20 seconds** before restarting the instrument in order to avoid damaging its electronics.

Chemical and biological safety

⚠ DANGER

Working with chemical samples, standards and reagents is linked with danger. It is recommended that the user of these products familiarize himself before the work with safe procedures and the correct use of the chemicals and to read all relevant safety data sheets carefully.

During the analysis of the samples, the use of toxic, highly flammable or corrosive chemicals can be required.

- Please observe safety information present on the packaging and material safety data sheets.
- All consumed solutions must be disposed of in accordance with the local and national regulations and laws.

Introduction

The additional Brewery Analysis software LZV936 is a compilation of all spectrophotometric applications that are relevant for brewery analysis. The methods are taken from the MEBAK user manuals and refer to the 4th edition published in 2002. For many analyses, the SIP 10 sipper offers the comfortable execution of the test with a pour-through cuvette.

Selection of a stored test

1. Select **STORED PROGRAMS**. An alphabetically arranged list of all available programs is shown. The methods of the additional Brewery Analysis software LZV936 are located at the end of the list.
2. Select a program, by choosing the corresponding test.
Note: *The programs of the Brewery Analysis have the numbers 2001 to 2017. Select **SELECT BY NUMBER** and confirm the entry with **OK**.*
3. Press **START** to start the program.

List of abbreviations

General notes

Reagents, if not otherwise noted, are of analytical purity.
Solutions, if not otherwise noted, are aqueous.

dist. H₂O	Distilled or demineralized water
sec	Seconds
min	Minutes
h	Hours
s	Standard deviation
r	Repeatability (95 %)
R	Comparability (95 %)
V_k	Variation coefficient
m	Mean value

Literature

MEBAK brew-technical analysis methods

Methods collection of the Mitteleuropäischen Brautechnischen Analysenkommission (MEBAK, Central European brew-technical analysis commission)

Published by the Chairman Dr. Heinrich Pfenninger

Self-publication of the MEBAK

D-85350 Freising-Weihenstephan

4th edition, 2002

Work procedures

Bitter units (EBC method)

Principle

Bitter substances, mainly iso- α acids, are extracted with iso-octane from the acidified sample and the concentration in the extract is determined spectrophotometrically.

Application area

Beer, worts

Measurement range

10–40 bitter units (beer)

20–60 bitter units (wort)

Accessories

- Centrifuge glasses with solvent-proof seals, 35 mL
- Glass beads
- Shaking machine
- Centrifuge, 3000 rpm
- Spectrophotometer, 275 nm
- 10 mm rectangular cuvette QS

Reagents

- Hydrochloric acid, 6N
- Iso-octane (2,2,4-trimethylpentane), spectroscopic pure (extinction measured in 10 mm rectangular cuvette QS at **275 nm** against dist. H₂O <0.010) (e.g. Uvasol)

Sample preparation

1. Separate wort and cloudy beer through centrifuging at **3000 rpm** for **15 min** (do not filter sample)
2. Remove carbon dioxide from the beer without removing any foam.
3. Pipette **10 mL** (or 5 mL wort + 5 mL dist. H₂O) of the sample maintained at **20 °C** (68 °F) into a centrifuge glass.
4. Add **0.5 mL** 6N hydrochloric acid, **20 mL** iso-octane and **three** glass beads.
5. Seal the centrifuge glass and leave to shake mechanically for **15 min** at **20 °C** (68 °F).
6. Centrifuge for **3 min** at **3000 rpm**.
7. Measure the extinction of the iso-octane extract in 10 mm rectangular cuvette at **275 nm** against iso-octane of the same quality (blank value).

Results specifications

Bitter units (BU) without decimals

Accuracy

$r = -0.36 + 0.05 m$

$R = 0.72 + 0.14 m$

Standard values

Beer: 10–40 BU, depending on nature, sort, type and origin

Wort: 20–60 BU, depending on beer and bitters yield

Literature

MEBAK brew-technical analysis methods
4th edition 2002, Volume II, P. 114 et seq

Procedure for bitter units

1. Prepare the blank value and the samples according to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2001** for the program "Beer bitter units". Confirm with **OK**.
Enter the number **2003** for the program "Wort bitter units". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the blank value (see sample preparation) and press **ZERO**.
6. Insert the analysis cuvette with the prepared iso-octane and press **READ**.

The result is shown on the display.

Note: In order to analyze further samples, repeat the procedure from step 6.

Total polyphenols (EBC method)

Principle

Polyphenols react with iron(III) ions in alkali solution under formation of colored iron complexes; the ensuing brownish coloring is measured spectrophotometrically.

Application area

Beer, worts

Measurement range

0–800 mg/L

Accessories

- Centrifuge
- Spectrophotometer, 600 nm
- 10 mm rectangular cuvette OS

Reagents

- Carboxymethylcellulose-ethylenediaminetetraacetic acid solution (CMC-EDTA-Na):
 - a. **Weigh out 10 g** CMC (low-viscosity) and 2 g EDTA-Na₂
 - b. Dissolve the substances by stirring in approximately **500 mL** H₂O, after complete dissolution fill with H₂O to **1 L**, if necessary, clear through centrifuging.
- Ammonium iron(III) citrate, 3.5%:
 - a. **Dissolve 3.5 g** ammonium iron(III) citrate, green (16 % Fe) in H₂O to **100 mL**; the solution must be completely clear and can be stored for approximately 1 week
 - b. Ammonia, diluted: dilute **1 part** concentrated ammonia (d = 0.91) with **2 parts** H₂O

Sample preparation

1. Remove carbon dioxide from the beer by shaking.
2. Clear the cloudy wort or beers through centrifuging.
3. Measure out **10 mL** sample + **8 mL** CMC-EDTA solution into a 25 mL volumetric flask and mix thoroughly.
4. Add **0.5 mL** iron(III) solution and mix **thoroughly**.
5. Add **0.5 mL** diluted ammonia solution and mix **thoroughly**.
6. Fill with H₂O to **25 mL** and mix **thoroughly**.
7. Measure the extinction after **10 min** in a 10 mm rectangular cuvette at **600 nm** against a blank value.
8. **Production of the blank value**
 - a. Measure out **10 mL** of sample (remove carbon dioxide in the beer with shaking, remove cloudiness of wort or beer with centrifuging) into a 25 mL volumetric flask.
 - b. Add **8 mL** CMC-EDTA solution and mix **thoroughly**.
 - c. Add **0.5 mL** diluted ammonia solution and mix **thoroughly**.
 - d. Fill with H₂O to **25 mL** and mix **thoroughly**.

Note: Mix the solution thoroughly after the addition of each individual reagent..

Results specifications

In mg/L without decimals

Accuracy

$s = \pm 9$

Standard values

Beer: 150–200 mg/L

Literature

MEBAK brew-technical analysis methods
4th edition 2002, Volume II, P. 107 et seq

Procedure for total polyphenols

1. Prepare the blank value and the samples according to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2002** for the program "Total polyphenols". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the blank value (see sample preparation) and press **ZERO**.
6. Insert the analysis cuvette with the prepared sample and press **READ**.

The result is shown on the display.

Note: *In order to analyze further samples, repeat the procedure from step 6.*

Reducibility (spectrophotometric method)

Reducibility is a measure for the fast-reducing substances present in the beer. Reducers occur in the beer in relatively small quantities, but have great significance in the chemical, physical, and biological properties as well as the taste stability of the beer.

Principle

Reducers reduce a certain amount of Tillmann's reagent (2,6-dichlorophenol indophenol, DPI) within a certain time span. The decoloration of the reagent is spectrophotometrically measured and calculated.

Measurement range

0–100

Accessories

- Spectrophotometer, 520 nm
- 10 mm rectangular cuvette OS
- Stop watch
- Water jet pump

Reagents

- 2,6-dichlorophenol-indophenol, 0.005 M (DPI solution, molecular weight of the sodium salt 290.08):
 - a. Weigh into a beaker approximately **100 mg** DPI, add approximately **25 mL** dist. H₂O and dissolve by heating (to approximately **60 °C** (140 °F)).
 - b. Rinse the contents after cooling into a 50 mL volumetric flask and fill to **50 mL**. Filter through a fiberglass filter.
 - c. Add **10 mL** filtrate, **1 g** KJ and **2 mL** H₂SO₄ (1+6) to a 150 mL Erlenmeyer flask and titrate with 0.01 N sodium thiosulfate until the color changes toward starch paste.
 - d. Calculate the indicator content:
used mL × 14.5 = mg indicator in 100 mL
 - e. Dilute the remaining filtrate so that 100 mL solution contains 145 mg indicator.
 - f. The solution can be stored in brim-full brown bottles at **+4 °C** (39 °F) for approximately 1 week.
- Phosphate-citrate buffer, pH =4.35:
- Dissolve **31.60 g** di-sodium hydrogen phosphate (Na₂HPO₄ × 12 H₂O) and **11.75 g** citric acid (C₆H₈O₇ × H₂O) in dist. H₂O and dilute to **1 L**.

Sample preparation

1. Maintain the beer at **20 °C** (68 °F) and remove the carbon dioxide under vacuum (water jet pump)
2. Pipette the decarbonated beer into a **10 mL** reagent flask with glass stopper. **Pour** the **0.25 mL** 0.005 M DPI solution slowly into the flask.
3. Seal the reagent flask **immediately** and invert the flask **two times**. Start the stop watch after the first inversion. After the first tip start the stop watch.
4. Fill a 10 mm rectangular cuvette with the mixture **immediately**. Measure the extinction against the blank value (decarbonated beer without the added reagent) at **520 nm** 60 seconds after the addition of the reagent.

Results specifications

Proportion of the supplied DPI quantity in %, which is reduced by 10 mL beer in 60 sec.

Accuracy

$$V_{kr} = \pm 1 \%$$

Standard values

> 60	very good
50–60	good
45–50	satisfactory
< 45	poor

Literature

MEBAK brew-technical analysis methods
4th edition 2002, Volume II, P. 104 et seq

Procedure for reducibility

1. Prepare the blank values and samples according to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2004** for the program "Reducibility". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the blank value (see sample preparation) and press **ZERO**.
6. Insert the analysis cuvette with the prepared sample. After 60 sec press **READ**.

The result is shown on the display.

Note: To analyze further samples, repeat the procedure from step 5, that means each sample requires a specific blank value.

Anthocyanogens (Harris and Ricketts method)

Anthocyanogens (leucoanthocyanidins) are phenolic compounds that are transformed by hot hydrochloric acid into red-colored anthocyanidins. Amount and condensation- or polymerization degree of these compounds have an influence on the formation of colloidal turbidities in the beer. Stabilization measures with PVPP correlate to a reduction of the anthocyanogen content.

Principle

The anthocyanogens are adsorbed to polyamide, the adsorbate dissolved in butanol-hydrochloric acid and heated, whereby a red-colored solution ensues, whose intensity is measured photometrically.

Application area

Beer, wort

Measurement range

0–100 mg/L

Accessories

- Shaking machine
- Centrifuge
- Mixing flask with ground-in stopper, 50 mL
- Frit 1 G4
- Filter flask
- Reagent flasks with ground-in stopper, 30 mL, graduation up to 25 mL
- Vacuum pump
- Spectrophotometer, 550 nm
- 10 mm rectangular cuvette OS

Reagents

- MN-polyamide SC 6
- **Solution 1:** n-Butanol / 37 % hydrochloric acid 5+1 (V/V)
- **Solution 2:** Dissolve **120 mg** iron(II) sulfate ($\text{FeSO}_4 \times 7 \text{H}_2\text{O}$) in 100 mL solution 1

Sample preparation

1. Centrifuge worts and young beers for **10 min** at **3000 rpm**.
2. Pipette **5 mL** beer (or wort) + **5 mL** dist. H_2O in a 50 mL mixing flask.
3. Pipette **10 mL** dist. H_2O (blank value) in a 50 mL mixing flask.
4. Rinse **0.5 g** polyamide powder with **10 mL** dist. H_2O in the mixing flask.
5. Mechanically shake the two mixing flasks for **40 min**.
6. Filter each of the suspensions through a 1 G4 frit and rinse **twice** with approximately **20 mL** dist. H_2O .
7. Vacuum the frits and polyamide powder dry, transfer the residue respectively with a spatula quantitatively into the reagent flask and purge with **15 mL** of solution 1.

8. Add **0.5 mL** of solution 2 to each flask and heat both reagent flasks for **30 min** in a boiling water bath (during the first 5 min stir well with a glass rod).
9. Take out the glass rods, rinse with a little solution 1, maintain the reagent flasks at 20 °C (68 °F) and fill with solution 1 to 25 mL.
10. Measure the extinction of the solution in a 10 mm rectangular cuvette at **550 nm** against an equivalently treated blank value (10 mL dist. H₂O instead of beer).

Results specifications

In mg/L without decimals

Accuracy

r = 9

Standard values

50–70 mg/L depending on raw materials and technological measures; with stabilization with PVPP correspondingly less.

Literature

MEBAK brew-technical analysis methods, 4th edition 2002, Volume II, P. 109 et seq

Procedure for anthocyanogens

1. Prepare the blank value and the samples according to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2005** for the program "Anthocyanogens". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the blank value (see sample preparation) and press **ZERO**.
6. Insert the analysis cuvette with the prepared sample and press **READ**.

The result is shown on the display.

Note: In order to analyze further samples, repeat the procedure from step 6.

Beer color (spectrophotometric EBC method)

Principle

This method eliminates the subjective influences of the human eye and differences in color impression with comparison of samples with a color standard. This technical method is considered the **official reference method**.

The extinction is measured in a 10 mm rectangular cuvette at a wavelength of 430 nm. The color in EBC units is derived by converting with a suitable factor.

Application area

Operating worts, beer, liquid malt substitutes of any type

Measurement range

0–60 units

Accessories

- Spectrophotometer, 430 nm (± 0.5 nm)
- 10 mm rectangular cuvette OS

Sample preparation

1. Dilute the sample so that the extinction lies within the linearity of the spectrophotometer.
2. Filter the sample through a membrane filter. The filtration can be omitted if the turbidity of the diluted sample is below 1 EBC turbidity unit.
3. Clarify the beer as required by addition of 0.1 % diatomaceous earth and a filter upstream of the membrane filter.
4. Measure the extinction (E) at **430 nm** against dist. H₂O (blank value).

Results specifications

In EBC units with 2 indicative numerals

Interferences

A spectrometric absorption curve does not reflect the color impression of the human eye, as light of equal intensity in various parts of the spectrum influences the eye differently. Also the extinction curves are at 430 nm very steep, so that slight measuring errors can occur. There are also differences during comparison of light beers with diluted dark beers.

Literature

MEBAK brew-technical analysis methods, 4th edition 2002, Volume II, P. 88 et seq

Procedure for beer color

1. Prepare the zero solution and the samples accordance to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2006** for the program "Beer color". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the zero solution (dist. H₂O) and press **ZERO**.
6. Insert the analysis cuvette with the prepared sample and press **READ**.

The result is shown on the display.

Note: In order to analyze further samples, repeat the procedure from step 6.

Free amino nitrogen (ninhydrin method according to EBC)

Principle

The examination solution is heated at pH 6.7 with ninhydrin and the resultant color measured at 570 nm. The method records the amino acids, ammonia and also the terminal alpha amino groups of peptides and proteins. Prolin is partially co-determined at the applied wavelength. The method is not specific to alpha-amino-nitrogen, because gamma-amino-butyric acid, which occurs in worts, also develops a color with ninhydrin.

Application areas

Beer, worts

Measurement range

0–400 mg/L

Accessories

- Reagent flasks with ground-in stoppers, 16 × 150 mm
- Variable pipette 1.0–5.0 mL (BBP065)
- Pipette tips for pipettes (BBP068)
- Boiling water bath
- Water bath at 20 °C (68 °F)
- Spectrophotometer, 570 nm
- 10 mm rectangular cuvette OS

Reagents

- **Color reagent:** Dissolve **10.0 g** di-sodium-hydrogen phosphate ($\text{Na}_2\text{HPO}_4 \times 12 \text{ H}_2\text{O}$), **6.0 g** potassium di-hydrogen phosphate (KH_2PO_4), **0.5 g** ninhydrin, **0.3 g** fructose in dist. H_2O and fill up to 100 mL (This solution can be stored for 2 weeks in a dark flask. The pH value must be 6.6–6.8).
- **Thinner:** Dissolve **2 g** potassium iodate in **600 mL** dist. H_2O and add **400 mL** 96 % ethanol
- **Parent solution:** Dissolve **107.2 mg** glycine in 100 mL dist. H_2O . Store this parent solution at 0 °C (32 °F).
- **Standard solution:** Fill up **1 mL** parent solution with dist. H_2O to **100 mL**. This standard solution contains 2 mg/L amino nitrogen.

Sample preparation

1. Dilute wort 100 times, beer 50 times (1–3 mg/L amino nitrogen).
2. Analyze the sample, standard solution and the blank value three times each.
3. Pipette **2 mL** of the diluted sample or the standard solution or the dist. H_2O into a reagent flask.
4. Add **1 mL** color reagent and mix.
5. Loosely seal the reagent flasks with glass stoppers to prevent loss due to evaporation.
6. Heat the reagent flasks for exactly **16 min** in a boiling water bath, then allow to cool for **20 min** in the water bath of 20 °C (68 °F).

7. Add **5 mL** of thinner.
8. Measure the extinction within 30 min in a 10 mm rectangular cuvette at **570 nm** against the equivalently treated blank value (dist. H_2O + color reagent).
9. **Correction for dark worts and beers (carry out three times)**
 - a. Add **2 mL** of the diluted sample to a reagent flask.
 - b. Add **1 mL** dist. H_2O instead of the color reagent and proceed as described above.
 - c. After adding **5 mL** of thinner measure against dist. H_2O .

Results specifications

In mg/L without decimals

Accuracy

$r = 17$
 $R = 28$

Standard values

Cast wort (12 %): 200–250 mg/L

Beer (12 %): 100–120 mg/L

For a satisfactory primary and secondary fermentation approximately 220–250 mg/L of free amino nitrogen should be present in the engaged wort.

Interferences

Because this test measures small amounts of amino acids, contamination must be avoided. The carefully cleaned flasks must only be touched on the outside. Ground-in stoppers etc. must be handled with tweezers.

Literature

MEBAK brew-technical analysis methods, 4th edition 2002
Volume II, P. 62 et seq

Comment

The following procedures describe a triplicate determination of blank values, standard solution, and samples without correction of light beer and worts.

For dark beers and worts, the work procedure contains a triplicate determination of blank value, standard solution, correction and sample.

Procedure for free amino nitrogen (FAN) in light worts and light beers

1. Prepare the zero solution, blank values, standard solutions and sample each three times according to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2007** for the program "FAN light wort". Confirm with **OK**.
Enter the number **2008** for the program "FAN light beer". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the zero solution (dist. H₂O) and press **ZERO**.
Z1 is shown on the display.
6. Insert the blank value cuvette (see sample preparation) and press **READ**.
R1 is shown on the display.
*Note: Proceed similarly with the blank value cuvettes 2 and 3. **R2** or **R3** are shown on the display.*
7. Insert the standard cuvette (see sample preparation) and press **READ**.
R4 is shown on the display.
*Note: Proceed similarly with the standard cuvettes 2 and 3. **R5** or **R6** are shown on the display.*

8. Insert the analysis cuvette with the prepared sample and press **READ**.

R7 is shown on the display.

*Note: Proceed similarly with the analysis cuvettes 2 and 3. **R8** or the result will be shown on the display.*

The FAN result is shown in mg/L.

Note: In order to analyze further samples, repeat the procedure from step 8.

Procedure for free amino nitrogen (FAN) in dark worts and dark beers

1. Prepare the zero solution, blank values, standard solutions, correction solutions and sample each three times according to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2015** for the program "FAN dark wort". Confirm with **OK**.
Enter the number **2016** for the program "FAN dark beer". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the zero solution (dist. H₂O) and press **ZERO**.
Z1 is shown on the display.
6. Insert the blank value cuvette (see sample preparation) and press **READ**.
R1 is shown on the display.
*Note: Proceed similarly with the blank value cuvettes 2 and 3. **R2** or **R3** are shown on the display.*
7. Insert the standard cuvette (see sample preparation) and press **READ**.
R4 is shown on the display.
*Note: Proceed similarly with the standard cuvettes 2 and 3. **R5** or **R6** are shown on the display.*

8. Insert the correction cuvette (see sample preparation) and press **READ**.
R7 is shown on the display.

Note: Proceed similarly with the correction cuvettes 2 and 3. **R8** or **R9** will be shown on the display.

9. Insert the analysis cuvette with the prepared sample and press **READ**.

R10 is shown on the display.

Note: Proceed similarly with the analysis cuvettes 2 and 3. **R11** or the result will be shown on the display.

The FAN result is shown in mg/L.

Note: In order to analyze further samples, repeat the procedure from step 9.

Steam-volatile phenols

The degree of smoking of whiskey malts is determined through the measurement of steam-volatile phenols. For beer production smoke malts are required only in small quantities for smoked beers, a Franconian specialty. Technical malfunctions during the drying process can however cause the malts, which for normal beers are used for processing, to demonstrate a smoky taste. The smoky taste can be detected by customers in beers that should not have that taste.

Next to the organoleptic test, the spectrophotometric determination of the steam-volatile phenols has proved the best method to establish which malt delivery produces the undesirable smoke taste, and to what extent this effects tank or bottled beers.

Principle

The phenol fraction extracted through steam is transformed with 4-amino-2,3-dimethyl-1-phenyl-3-pyrazolin-5-one (4-aminophenazone) in the alkali environment and under the oxidation effect of potassium hexacyanoferrate(III) into a pigment, which can be measured after chloroform extraction spectrophotometrically.

Application areas

Malt, beer

Measurement range

0–20 mg/kg

Comments

Wheat beers cannot be analyzed using this method because, due to the activity of the top-fermented yeast, they contain a high quantity of steam-volatile phenols but no smoke flavor.

Accessories

- DLFU mill, mill gap 1 mm
- Steam distillation unit
- Separatory funnel, 1 L
- Spectrophotometer, 460 nm
- 40 mm rectangular cuvette OS

Reagents

- Chloroform, high purity
- Silicone anti-foam emulsion
- Phosphoric acid, concentrated (d=1.71)
- Copper sulfate, $\text{CuSO}_4 \times 5 \text{H}_2\text{O}$, 10 %
- Ammonium chloride, 5 %
- 4-amino-2,3-dimethyl-1-phenyl-3-pyrazolin-5-one, 2 %: prepare fresh daily
- Potassium hexacyanoferrate(III), $\text{K}_3[\text{Fe}(\text{CN})_6]$, 8 %: prepare fresh daily

- Phenol standard solution: Dissolve 1.000 g phenol in dist. H₂O to 1000 mL (1 mL=1 mg). The solution must be clear and colorless. From this solution freshly prepare thinners for constructing the calibration curve 0.02–0.1 mg/L as required.
- Dilute ammonia, diluted (1+4): 1 part concentrated ammonia (d=0.91) with 4 parts dist. H₂O.

Sample preparation

1. Steam distillation

- Add **50 g** coarse malt with **500 mL** dist. H₂O (for tests of 300 mL of beer) to a distillation flask.
- Add **3 mL** copper sulfate solution.
- Adjust the pH value through addition of phosphoric acid to pH<4.
- Add the silicone anti-foam emulsion.
- Perform steam distillation until 300 mL have been recovered.

2. Color reaction

- Add **10 mL** of Ammonium Chloride solution to the distillate. For real smoke malts or whiskey malts, use correspondingly less distillate (for example 100 mL).
To produce the blank value, introduce 300 mL dist. H₂O instead of distillate and also add **10 mL** ammonium chloride solution.
- Shake the solution.
- In each case (distillate and blank value) adjust the pH value with ammonia to pH 10.2 ± 0.1.
- Transfer each solution to a 1 L separatory funnel.
- Add to each **3 mL** 4-amino-2,3-dimethyl-1-phenyl-3-pyrazolin-5-one and **3 mL** potassium hexacyanoferrate(III).
- Shake the solution.
- Allow the solutions to rest for **3 min**.
- Shake each separatory funnel **3 times at one minute each** with **10 mL** chloroform to extract.
- Wait **10 min** each for the phase separation.

- Filter the chloroform extracts through paper filters into 25 mL volumetric flasks.
- Rinse the filters with some chloroform.
- Fill with chloroform respectively up to the mark.
- Measure the chloroform extract (prepared distillate) in 40 mm rectangular cuvette at **460 nm** against the prepared blank value.

3. Calibration values

- With phenol standard solutions of 0.02–0.1 mg/L (introduce 300 mL) carry out steam distillation and convert as described above.

Results specifications

In mg/kg with two decimals (or mg/L for beer)

Accuracy

Vk = ± 5 % (repeat error)

Target values

Malts: < 0.2 mg/kg: no smoke taste expected

Beers: < 0.03 mg/L: harmless in most cases

The penetration of the smoke taste depends on the beer composition, which is why the named lower limit applies only with reservations.

Literature

MEBAK brew-technical analysis methods, 3rd edition, Volume I

Procedure for steam-volatile phenols

- Prepare the blank value and the samples according to the work procedure.
- Select **STORED PROGRAMS > SELECT BY NUMBER**.
- Enter the number **2009** for the program "Steam-volatile phenols". Confirm with **OK**.
- Press **START** to start the program.
- Enter the blank value (see sample preparation) and press **ZERO**.

6. Insert the analysis cuvette with the prepared sample and press **READ**.

The result is shown on the display.

Note: *In order to analyze further samples, repeat the procedure from step 6.*

Photometric iodine sample

Photometric iodine sample according to the new method (MEBAK from 1993)

Principle

High-molecular dextrans and starches are precipitated through addition of ethanol to wort and beer, centrifuged out, dissolved in phosphate buffer and displaced with iodine solution. Depending on molecular weight and branching factor of the erythrodextrin and starch, a red to blue color forms, whose intensity is measured photometrically.

Application area

Wort, beer (samples with an iodine value > 0.8 must be diluted.)

0–1 iodine value

Accessories

- Centrifuge
- Centrifuge flasks with ground-in stoppers, 100–110 mL capacity
- Shaking machine
- Pipettes, 0.5 mL, 2 mL, 10 mL, 20 mL, 40 mL
- Spectrophotometer, 578 nm
- 40 mm rectangular cuvette OS
- "Plumper" or plastic spatula
- Reagents
- Ethanol, 95 %
- Iodine solution, 1N (parent solution)
- Iodine solution, 0.02N (prepare fresh daily from the parent solution)

- Phosphate buffer solution, 0.1M, pH 3.5: adjust 0.1M KH_2PO_4 solution with 0.1M phosphoric acid to pH 3.5

Sample preparation

1. Pipette **10.0 mL** centrifuged wort or decarbonated beer into a centrifuge flask.
2. Add **40.0 mL** ethanol and shake mechanically for **10 min**.
3. Centrifuge for **5 min** at **2500 rpm**.
4. Decant the clear phase as carefully and completely as possible.
5. Loosen the residue with **20.0 mL** phosphate buffer solution by mechanically shaking for **10 min**.
6. Centrifuge the solution for **5 min** at **2500 rpm**.
7. Pipette **2 mL** of the excess and **8 mL** phosphate buffer solution into a 40 mm rectangular cuvette and measure at **578 nm** against phosphate buffer solution.
8. Add **0.5 mL** 0.02N iodine solution, **immediately** mix the contents with the "spatula", measure after **30 sec**.
9. **Iodine blank value**
 - a. Pipette **10 mL** phosphate buffer solution and **0.5 mL** 0.02N iodine solution into a 40 mm rectangular cuvette and mix.
 - b. Measure the extinction at **578 nm** against phosphate buffer solution.

Results specifications

Extinction to 2 decimals

Accuracy

$V_{kr} = \pm 3 \%$

Standard values

< 0.3 (wort)

Literature

MEBAK brew-technical analysis methods, 4th edition 2002, Volume II, P. 34 et seq

Procedure for photometric iodine sample

1. Prepare the blank value (phosphate buffer), the iodine blank value and the samples in accordance with the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2010** for the program "Photometric iodine sample". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the blank value (phosphate buffer) and press **ZERO**.
Z1 is shown on the display.
6. Insert the iodine blank value cuvette (see sample preparation) and press **READ**.
R1 is shown on the display.
7. Insert the analysis cuvette with the prepared sample and press **READ**.
R2 is shown on the display.
8. Remove the analysis cuvette and dispense **0.5 mL** 0.02N iodine solution into the analysis cuvette.
9. Mix the content **immediately** with the "spatula" and insert the analysis cuvette after 30 sec. Press **READ**.
The result is shown on the display.

Note: In order to analyze further samples, repeat the procedure from step 7. The prepared iodine blank value solution can be used for the complete measurement series.

Thiobarbituric acid number (TAN)

The thiobarbituric acid number applies as the summary variable for the thermal load of malt and wort. It is an indicator that in addition to 5-hydroxymethylfurfural (HMF) records a multitude of products of the Maillard reaction and other organic compounds.

Principle

The examination sample is used with acetic thiobarbituric acid solution in reaction and the resulting yellow coloration is measured spectrophotometrically.

Application area

Beer, wort, congress wort or malt extract

Measurement range

0–100

Accessories

- Water bath, 70 °C (158 °F)
- Brown reaction flasks with ground-in stopper, 20 mL or 25 mL
- Spectrophotometer, 448 nm
- 10 mm rectangular cuvette OS

Reagents

- Acetic acid, 90 %:
dilute 225 g acetic acid 100 % (glacial acetic acid) with dist. H₂O to 250 g.
- Thiobarbituric acid, 0.02 mol/L: Dissolve **0.288 g** 2-thiobarbituric acid (M = 144.15 g/mol) in 100 mL volumetric flask with 90 % acetic acid by warming in the water bath and after cooling to 20 °C (68 °F) fill up to the mark with 90 % acetic acid. Prepare fresh daily.
- Diatomaceous Earth

Sample preparation

Note: Follow procedure carefully as accuracy can affect final results.

1. Clear cloudy examination solutions through filtration over diatomaceous earth.

2. Dilution

- Dilute worts and beers 10 times with dist. H₂O.
- Dilute congress worts 5 times with dist. H₂O.

3. Empty value

- a. **Displace 10 mL** of the diluted sample with **5 mL** 90 % acetic acid, shake and proceed as for the main value.

4. Main value

- a. **Displace 10 mL** of the diluted sample with **5 mL** thiobarbituric acid and shake.
- b. **Place** in a water bath of 70 °C for 70 min (158 °F) (avoid direct sunlight, make sure that when positioning the reagent flasks the temperature in the bath only temporarily sinks by 1–2 °).
- c. After elapse of the reaction time quickly cool down reagent flasks to 20 °C (68 °F) (fast-flowing cold water or cold bath).
- d. Measure yellow coloration **immediately** in a 10 mm rectangular cuvette at **448 nm** against dist. H₂O.

Results specifications

Thiobarbituric acid number (TAN), dimensionless number

Standard values

Light cast wort < 45

Light cold wort (after wort cooling) < 60

Literature

MEBAK brew-technical analysis methods, 4th edition 2002, Volume II, P. 35 et seq

Procedure for thiobarbituric acid number in beer/wort and congress wort

1. Prepare the zero solution, empty values and samples in accordance with the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2011** for the program "TAN C-wort". Confirm with **OK**.
Enter the number **2012** for the program "TAN beer/wort". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the zero solution (dist. H₂O) and press **ZERO**.
Z1 is shown on the display.
6. Insert the empty value cuvette (see sample preparation) and press **READ**.
R1 is shown on the display.
7. Insert the analysis cuvette with the prepared sample and press **READ**.

The result is shown on the display.

Note: In order to analyze further samples, repeat the procedure from step 6.

Iso- α and α acids

Principle

Bitters are extracted with iso-octane from the acidified sample and certain disruptive substances removed through washing of the extract with acidic methanol. The concentration of iso- α acids and α acids is determined through measuring of the extinction in alkali methanol at 255 nm and 360 nm.

Application area

Beer, wort

Measurement range

0–60 mg/L

Accessories

- Centrifuge flasks with solvent-proof screw closure, 100–110 mL capacity
- Shaking machine
- Centrifuge, 3000 rpm
- Spectrophotometer, 255 nm and 360 nm
- 10 mm rectangular cuvette QS

Reagents

- Hydrochloric acid, 6N
- Iso-octane (2,2,4-trimethylpentane), spectroscopic purity
- Sodium sulfate, anhydrous
- Methanol
- Hydrochloric acid, 4N
- Sodium hydroxide, 6N, decarbonated
- Acidic methanol solution: mix 64 mL methanol and 36 mL 4N hydrochloric acid (prepare fresh daily)
- Alkali methanol solution: top up 0.2 mL 6N sodium hydroxide with methanol to 100 mL (prepare fresh daily)

Sample preparation

1. Clear wort or cloudy beer through centrifuging at **3000 rpm** for **15 min** (do not filter).
2. Remove the carbon dioxide from the beer without foam loss.
3. Pipette **50 mL** of the sample maintained at 20 °C (68 °F) into a centrifuge flask.
4. Add **3 mL** 6N hydrochloric acid and **25 mL** iso-octane.
5. Seal the centrifuge flask and shake it mechanically for **30 min** at optimal mixing intensity.
6. Centrifuge to separate the phases and break up the emulsion for **5 min** at **3000 rpm**.

7. Remove the lower aqueous phase through suction with a pipette and discard it. Displace the iso-octane phase with enough sodium sulfate so that the solution is clear after brief vigorous shaking.
8. Pipette **10 mL** of the iso-octane phase into a 25 mL volumetric flask.
9. Add **10 mL** acidic methanol solution, seal the flask and mix the contents by inverting the flask 100 times.
10. Pipette **5 mL** of the excess clear iso-octane phase into a 25 mL volumetric flask.
11. Fill to the mark with alkali methanol solution and mix.
12. Measure the extinction of the iso-octane solution against a blank value at **255 nm** and **360 nm**.
13. **Blank value preparation**
 - a. Pipette **5 mL** iso-octane into a 25 mL volumetric flask.
 - b. Fill to the mark with alkali methanol solution and mix.

Results specifications

In mg/L without decimals

Accuracy

$V_{kr} = \pm 5 \%$

Standard values

Beer: 10–40 mg/L iso- α acids, depending on nature, sort, type and origin (< 2mg/L α acids)

Wort: 15–50 mg/L iso- α acids, depending on beer and bitters yield
1–15 mg/L α acids depending on degree of isomerization

Literature

Brew-technical analysis methods, 4th edition 2002, Volume II, P. 116 et seq

Procedure for iso- α and α acids

1. Prepare the blank value and the samples according to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2013** for the program "Iso- α and α acids". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the blank value (see sample preparation) and press **ZERO**.
6. Insert the analysis cuvette with the prepared sample and press **READ**.

The result is shown on the display.

Note: In order to analyze further samples, repeat the procedure from step 6. The prepared blank value can be used for the complete measuring series.

Vicinal diketones (diacetyl, 2,3-pentanedione)

During yeast metabolism, fermentation causes 2-acetolactate and 2-acetohydroxy butyrate to emerge. Oxidation causes a conversion into the vicinal diketones diacetyl (2,3 butanedione) and 2,3-pentanedione. Diacetyl can however also occur as a characteristic metabolic product of certain micro-organisms. As soon as the threshold value is exceeded, the beer obtains a bad aroma.

The photometric determination method in operation control is frequently preferred over the gas chromatographic methods, as it can be performed quickly and without great equipment expense. It does not however permit the preferable differentiation between diacetyl and pentanedione.

Principle

The basis of the method is the reaction between diacetyl and/or 2,3-pentanedione and 1,2-phenylenediamine under formation of 2,3-dimethylquinoxaline, which exhibits a specific absorption at 335 nm.

Application area

Beer

Measurement range

0–1 mg/kg

Accessories

- Appliance for nitrogen determination with heating jacket, macro execution (for example from Schott). The supplied cooler must possibly be replaced by a larger one if the distillate is not sufficiently cooled. Other similar appliances, for example from Büchi, are also suitable.
- Spectrophotometer, 335 nm
- 20 mm rectangular cuvette QS

Reagents

- Hydrochloric acid, 4N
- 1,2-phenylenediamine, 1 % in 4N hydrochloric acid (prepare solution fresh on the day of application and store in the dark). 1,2-phenylenediamine is toxic and an allergen; it must be handled carefully, work with gloves
- Anti-foam emulsion (free from diketones)

Sample preparation

1. Add **100 g** non-decarbonated beer into a pre-heated distillation appliance.
2. Add a drop of anti-foam emulsion.
3. Control the steam supply so that within 2 min approximately 25 mL distillate transfers.
4. Catch the distillate in a 25 mL volumetric flask.
5. Pipette **10 mL** of the mixed distillate into two 50 mL Erlenmeyer flasks (main value, blank value).
6. **Blank value**
 - Add **2.5 mL** 4N hydrochloric acid.

7. Main value

- Add **0.5 mL** 1,2-phenylenediamine solution, mix and place in the dark for **30 min**.
- Add **2 mL** 4N hydrochloric acid.

8. Measure the extinction of the main value against the blank value within **20 min** at **335 nm** in a 20 mm rectangular cuvette.

Results specifications

In mg/kg with two decimals

Accuracy

s = ± 0.01

Target value

For light full beer <0.15 mg/kg

Comments

The distillation appliance does not need to be cleaned or purged between individual tests and can be immediately refilled between tests. After testing is complete, the adherent residues should be removed with dilute sodium hydroxide or another suitable cleaning agent.

Acetohydroxy acids present in bottled beer are oxidized in the presence of oxygen to diketones. The beer sample can, for analysis of the total diketone content before the actual analysis, be maintained for 1.5 h at 70 °C (158 °F).

Literature

MEBAK brew-technical analysis methods, 4th edition 2002
Volume II, P. 134 et seq

Procedure for vicinal diketones

1. Prepare the blank value and the samples according to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2014** for the program "Vicinal diketones". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the blank value (see sample preparation) and press **ZERO**.
6. Insert the analysis cuvette with the prepared main value and press **READ**.

The result is shown on the display.

Note: To analyze further samples, repeat the procedure from step 5.

Iron

Iron can dissolve into beer through raw materials, filter additives, and clarifiers. It can also be absorbed from appliances, lines, containers, and beer foam stabilizing agents. Iron negatively effects the colloidal stability, the taste and the gushing tendency of the beer.

Principle

Bivalent iron forms a violet colored complex with a very high molar extinction coefficient with the di-sodium salt of 5,6-diphenyl-3-(2-pyridyl)-1,2,4-triazine-4,4-disulfonic acids (ferrozine). Trivalent iron must be reduced to bivalent iron before measurement. The color intensity is measured spectrometrically.

Measurement range

0–1 mg/L

Accessories

- Weighing balances, reading precision 0.1 mg
- Pipettes, 0.1 mL, 2 mL, 5 mL
- Spectrophotometer, 560 nm
- 40 mm rectangular cuvette OS

Reagents

Prepare all solutions with iron-free dist. H₂O.

- Buffer solution, pH 4.3: dissolve **75 g** ammonium acetate and **150 g** concentrated acetic acid in approx. **800 mL** dist. H₂O, check pH and top up to 1 L.
- Ferrozine reagent: dissolve **0.257 g** ferrozine or ferrospectral in **50 mL** buffer (solution can be stored for 2 weeks)
- Ascorbic acid, 2.5 % (prepare fresh daily)
- Hydrochloric acid, concentrated
- Iron(III) standard solution for determination of the calibration curve: dissolve **863.4 mg** ammonium iron(III) sulfate [NH₄Fe(SO₄)₂ × 12 H₂O] in the 1 L volumetric flask in dist. H₂O and after addition of **0.1 mL** concentrated hydraulic acid top up with dist. H₂O to the mark. After dilution of 50 mL of this solution with dist. H₂O to 1 L, the resulting standard solution contains 5 mg/mL Fe³⁺.

Sample preparation

1. Remove the carbon dioxide in the beer and allow the foam to subside completely.
2. Pipette **40 mL** beer, **2 mL** ferrozine reagent, and **1 mL** ascorbic acid solution into a 50 mL volumetric flask.
3. Fill to the mark with dist. H₂O.
4. Prepare the blank value without addition of ferrozine reagent. Use a separate blank for each individual beer.
5. Measure the extinction of the solution in a 40 mm rectangular cuvette at **560 nm** against a corresponding blank value.

Results specifications

In mg/L with three significant places

Accuracy

r = 0.0080

Target value

< 0.200 mg/L

Literature

MEBAK brew-technical analysis methods, 4th edition 2002
Volume II, P. 149 et seq

Procedure for iron

1. Prepare the blank values and samples according to the work procedure.
2. Select **STORED PROGRAMS > SELECT BY NUMBER**.
3. Enter the number **2017** for the program "Iron". Confirm with **OK**.
4. Press **START** to start the program.
5. Enter the blank value (see sample preparation) and press **ZERO**.
6. Insert the analysis cuvette with the prepared main value and press **READ**.

The result is shown on the display.

Note: To analyze further samples, repeat the procedure from step 5.

Determination of the calibration curve

The factor $1 = 0.037$ is an empirical factor and must be individually determined through a calibration curve. The factor is the slope of the calibration curves.

1. Pipette **40 mL** of beer into four individual 50 mL volumetric flasks.
2. Pipette into the volumetric flasks each **0.40 mL**, **0.80 mL**, **1.60 mL** and **3.20 mL** of the iron standard solution (5 mg Fe^{3+} /mL).
3. Add to each volumetric flask **2 mL** ferrozine reagent and **1 mL** ascorbic acid solution.
4. Fill to the mark with dist. H_2O .
5. Measure the extinction of the solution in a 40 mm rectangular cuvette at 560 nm against a corresponding blank value.

Subtract the results of the un-spiked sample from the spiked samples.

Using these programs with the SIP 10 sipper

Information on the installation, configuration, and sample insertion when using the SIP 10 sipper can be found in the *SIP 10 user manual*.

Accessories

Description	Order no.
Rectangular cuvettes set (standard deviation = 10 mm, aligned pair)	2095100
Rectangular cuvette QS, standard deviation = 10 mm (3.5 mL)	2624410
Rectangular cuvette OS, standard deviation = 10 mm (3 units)	LZP045
Rectangular cuvette QS, standard deviation = 20 mm	LZV008
Rectangular cuvette OS, standard deviation = 20 mm	LZP331
Rectangular cuvette QS, standard deviation = 50 mm (17.5 mL)	2624450
Rectangular cuvette OS, standard deviation = 50 mm	LZP167
Pour-through cuvette QS, Standard deviation = 3 mm, center height = 10 mm, height = 40 mm	LZV638
Pour-through cuvette QS, Standard deviation = 50 mm, center height = 10 mm, height = 40 mm	LZV649
Pour-through cuvette QS, Standard deviation = 10 mm, center height = 10 mm, height = 40 mm	LZV510
SIP 10 sipper module set for DR 6000 complete with tray, tube set and 1 cm quartz glass pour-through cuvette, EU	LQV157.99.30001

**HACH Company
World Headquarters**

P.O. Box 389
Loveland, Colorado
80539-0389 U.S.A.
Tel (800) 227-HACH
(800) -227-4224
(U.S.A. only)
Fax (970) 669-2932
orders@hach.com
www.hach.com

HACH LANGE GMBH

Willstätterstraße 11
D-40549 Düsseldorf
Tel. +49 (0)2 11 52 88-320
Fax +49 (0)2 11 52 88-210
info@hach-lange.de
www.hach-lange.de